




GFM2 Antibody

Product Code	CSB-PA924833
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q969S9
Immunogen	Synthesized peptide derived from internal of Human GFM2.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total GFM2 protein.
Tested Applications	ELISA,WB,IF;WB:1:500-1:3000,IF:1:100-1:500
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	EF-G 2; EFG2; elongation factor G 2; mitochondrial; elongation factor G2
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	GFM2

Image



Western blot analysis of extracts from LOVO cells and RAW264.7 cells, using GFM2 antibody.



Immunofluorescence analysis of A549 cells, using GFM2 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.