




ELOVL1 Antibody

Product Code	CSB-PA262320
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BW60
Immunogen	Synthesized peptide derived from internal of Human ELOVL1.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Specificity	The antibody detects endogenous levels of total ELOVL1 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Elongation of very long chain fatty acids protein 1; ELOVL1; SSC1; CGI-88;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ELOVL1

Image



Western blot analysis of extracts from HepG2 cells and Jurkat cells, using ELOVL1 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.