





ELOVL5 Antibody

Product Code	CSB-PA870663
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NYP7
Immunogen	Synthesized peptide derived from internal of Human ELOVL5.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total ELOVL5 protein.
Tested Applications	ELISA,WB,IHC;WB:1:500-1:3000,IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Elongation of very long chain fatty acids protein 5; ELOV5; ELOVL2; hELO1;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ELOVL5

Image



Western blot analysis of extracts from MCF-7 cells, using ELOVL5 antibody.



Immunohistochemistry analysis of paraffin-embedded human brain tissue using ELOVL5 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.