




NEIL3 Antibody

Product Code	CSB-PA065468
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TAT5
Immunogen	Synthesized peptide derived from internal of Human NEIL3.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total NEIL3 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Endonuclease VIII-like 3; Nei-like 3; DNA glycosylase FPG2;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	NEIL3

Image



Western blot analysis of extracts from COLO cells, HUVEC cells, HepG2 cells and HeLa cells, using NEIL3 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.