




MAN1B1 Antibody

Product Code	CSB-PA217897
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UKM7
Immunogen	Synthesized peptide derived from internal of Human MAN1B1.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Liquid
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MAN1B1

Image



Western blot analysis of extracts from Hela cells, using MAN1B1 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.