





EPN3 Antibody

Product Code	CSB-PA243224
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H201
Immunogen	Synthesized peptide derived from internal of Human EPN3.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total EPN3 protein.
Tested Applications	ELISA, WB, IHC; WB:1:500-1:3000, IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Epsin-3; EPS-15-interacting protein 3;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	EPN3

Image



Western blot analysis of extracts from K562 cells, using EPN3 antibody.



Immunohistochemistry analysis of paraffin-embedded human colon carcinoma tissue, using EPN3 antibody.