




EIF3K Antibody

Product Code	CSB-PA693741
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UBQ5
Immunogen	Synthesized peptide derived from internal of Human EIF3K.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total EIF3K protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Eukaryotic translation initiation factor 3 subunit K; eukaryotic translation initiation factor 3 subunit 12; eIF-3 p25; eIF-3 p28; eIF3k; Muscle-specific gene M9 protein
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	EIF3K

Image



Western blot analysis of extracts from HeLa cells, using EIF3K antibody.