




# EIF3K Antibody

<b>Product Code</b>	CSB-PA693741
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q9UBQ5
<b>Immunogen</b>	Synthesized peptide derived from internal of Human EIF3K.
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Specificity</b>	The antibody detects endogenous levels of total EIF3K protein.
<b>Tested Applications</b>	ELISA, WB; WB: 1:500-1:3000
<b>Form</b>	Rabbit IgG in phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
<b>Purification Method</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Eukaryotic translation initiation factor 3 subunit K; eukaryotic translation initiation factor 3 subunit 12; eIF-3 p25; eIF-3 p28; eIF3k; Muscle-specific gene M9 protein
<b>Product Type</b>	Polyclonal Antibody
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	EIF3K

## Image



Western blot analysis of extracts from HeLa cells, using EIF3K antibody.

**Usage** For Research Use Only. Not for use in diagnostic or therapeutic procedures.