




FKBPL Antibody

Product Code	CSB-PA904562
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UIM3
Immunogen	Synthesized peptide derived from internal of Human FKBPL.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total FKBPL protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	DIR1; FK506 binding protein like; WAF-1/CIP1 stabilizing 39; WISP39;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	FKBPL

Image



Western blot analysis of extracts from Jurkat cells, using FKBPL antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.