





GPRIN2 Antibody

Product Code	CSB-PA111091
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60269
Immunogen	Synthesized peptide derived from N-terminal of Human GPRIN2.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total GPRIN2 protein.
Tested Applications	ELISA, WB, IHC; WB: 1:500-1:3000, IHC: 1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	G protein-regulated inducer of neurite outgrowth 2; GRIN2;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	GPRIN2

Image



Western blot analysis of extracts from HuvEc cells, K562 cells and HepG2 cells, using GPRIN2 antibody.



Immunohistochemistry analysis of paraffin-embedded human brain tissue using GPRIN2 antibody.