




GABRA4 Antibody

Product Code	CSB-PA009683
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P48169
Immunogen	Synthesized peptide derived from internal of Human GABRA4.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total GABRA4 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Gamma-aminobutyric acid receptor subunit alpha-4; GABA(A) receptor subunit alpha-4;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	GABRA4

Image



Western blot analysis of extracts from HepG2 cells, using GABRA4 antibody.