




GGH Antibody

Product Code	CSB-PA789575
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92820
Immunogen	Synthesized peptide derived from internal of Human GGH.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total GGH protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	conjugase; EC 3.4.19.9; folylpolygammaglutamyl hydrolase; gamma-Glu-X carboxypeptidase; gamma-glutamyl hydrolase
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	GGH

Image



Western blot analysis of extracts from RAW264.7 cells, using GGH antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.