




IFI16 Antibody

Product Code	CSB-PA586413
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q16666
Immunogen	Synthesized peptide derived from internal of Human IFI16.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total IFI16 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Gamma-interferon-inducible protein I γ -16; Interferon-inducible myeloid differentiation transcriptional activator; IFI 16;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	IFI16

Image



Western blot analysis of extracts from HeLa cells and HepG2 cells, using IFI16 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.