




GJB7 Antibody

Product Code	CSB-PA449513
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6PEY0
Immunogen	Synthesized peptide derived from internal of Human GJB7.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total GJB7 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	connexin25; CX25; CXB7; gap junction protein; beta 7
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	GJB7

Image



Western blot analysis of extracts from HT-29 cells, using GJB7 antibody.