




GCFC2 Antibody

Product Code	CSB-PA578480
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P16383
Immunogen	Synthesized peptide derived from internal of Human GCF.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total GCF protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	C2orf3; GC-rich sequence DNA-binding factor; GCF; TCF-9; Transcription factor 9
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	GCFC2

Image



Western blot analysis of extracts from RAW264.7 cells, using GCF antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.