




TSTA3 Antibody

Product Code	CSB-PA790986
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13630
Immunogen	Synthesized peptide derived from internal of Human TSTA3.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total TSTA3 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	GDP-L-fucose synthetase; EC 1.1.1.271; Protein FX; Red cell NADP(H)-binding protein; GDP-4-keto-6-deoxy-D-mannose-3
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TSTA3

Image



Western blot analysis of extracts from rat kidney cells and rat muscle cells, using TSTA3 antibody.