




# TSTA3 Antibody

|                            |  |
|----------------------------|--|
| <b>Product Code</b>        | CSB-PA790986   |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.  |
| <b>Uniprot No.</b>         | Q13630   |
| <b>Immunogen</b>           | Synthesized peptide derived from internal of Human TSTA3.  |
| <b>Raised In</b>           | Rabbit   |
| <b>Species Reactivity</b>  | Human  |
| <b>Specificity</b>         | The antibody detects endogenous levels of total TSTA3 protein.   |
| <b>Tested Applications</b> | ELISA,WB;WB:1:500-1:3000   |
| <b>Form</b>                | Rabbit IgG in phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| <b>Purification Method</b> | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.                              |
| <b>Clonality</b>           | Polyclonal   |
| <b>Alias</b>               | GDP-L-fucose synthetase; EC 1.1.1.271; Protein FX; Red cell NADP(H)-binding protein; GDP-4-keto-6-deoxy-D-mannose-3                                |
| <b>Product Type</b>        | Polyclonal Antibody  |
| <b>Immunogen Species</b>   | Homo sapiens (Human)   |
| <b>Target Names</b>        | TSTA3  |

## Image



Western blot analysis of extracts from rat kidney cells and rat muscle cells, using TSTA3 antibody.

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.