




GSTT1/GSTT4 Antibody

Product Code	CSB-PA941427
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P30711
Immunogen	Synthesized peptide derived from internal of Human GSTT1/4.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total GSTT1/4 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	EC 2.5.1.18; glutathione S-transferase theta 1; glutathione transferase T1-1; GST class-theta-1;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	GSTT1

Image



Western blot analysis of extracts from Jurkat cells, using GSTT1/4 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.