




GLRB Antibody

Product Code	CSB-PA088748
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P48167
Immunogen	Synthesized peptide derived from internal of Human GLRB.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total GLRB protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Glycine receptor subunit beta; Glycine receptor 58 kDa subunit;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	GLRB

Image



Western blot analysis of extracts from A549 cells and K562 cells, using GLRB antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.