




GORASP2 Antibody

Product Code	CSB-PA446729
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H8Y8
Immunogen	Synthesized peptide derived from internal of Human GORASP2.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total GORASP2 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	golgi phosphoprotein 6; golgi reassembly-stacking protein 2; GOLPH6; GORS2; GRASP55
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	GORASP2

Image



Western blot analysis of extracts from HeLa cells, HepG2 cells, Jurkat cells and LOVO cells, using GORASP2 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.