




GRAP2 Antibody

Product Code	CSB-PA443192
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75791
Immunogen	Synthesized peptide derived from internal of Human GRAP2.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total GRAP2 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	GRB2-related adapter protein 2??ADS protein; Growth factor receptor-binding protein; GRBLG; Grf40 adapter protein; Grf-40
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	GRAP2

Image



Western blot analysis of extracts from HuvEc cells, using GRAP2 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.