




GADD45GIP1 Antibody

Product Code	CSB-PA039959
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TAE8
Immunogen	Synthesized peptide derived from internal of Human GADD45GIP1.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total GADD45GIP1 protein.
Tested Applications	ELISA,IHC;IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Growth arrest and DNA-damage-inducible proteins-interacting protein 1; CR6-interacting factor 1; CRIF1; CKII beta-associating protein; Papillomavirus L2-interacting nuclear protein 1
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	GADD45GIP1

Image



Immunohistochemistry analysis of paraffin-embedded human brain tissue using GADD45GIP1 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.