




DIRAS1 Antibody

Product Code	CSB-PA550322
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95057
Immunogen	Synthesized peptide derived from C-terminal of Human DIRA1.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total DIRA1 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	GTP-binding protein Di-Ras1; Distinct subgroup of the Ras family member 1; Ras-related inhibitor of cell growth; Rig protein; Small GTP-binding tumor suppressor 1
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	DIRAS1

Image



Western blot analysis of extracts from HeLa cells and HepG2 cells, using DIRA1 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.