




GNA14 Antibody

Product Code	CSB-PA998160
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95837
Immunogen	Synthesized peptide derived from N-terminal of Human GNA14.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total GNA14 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Guanine nucleotide-binding protein alpha-14 subunit; G-protein alpha subunit 14; GNA14;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	GNA14

Image



Western blot analysis of extracts from HuvEc cells, using GNA14 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.