




GNG5 Antibody

Product Code	CSB-PA032886
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P63218
Immunogen	Synthesized peptide derived from internal of Human GNG5.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Specificity	The antibody detects endogenous levels of total GNG5 protein.
Tested Applications	ELISA, IF; IF: 1:100-1:500
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Guanine nucleotide-binding protein G(I)/G(S)/G(O) subunit gamma-5;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	GNG5

Image



Immunofluorescence analysis of A549 cells, using GNG5 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.