





GNAL Antibody

Product Code	CSB-PA111285
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P38405
Immunogen	Synthesized peptide derived from internal of Human GNAL.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total GNAL protein.
Tested Applications	ELISA,WB,IHC;WB:1:500-1:3000,IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	adenylate cyclase-stimulating G alpha protein; olfactory type; guanine nucleotide-binding protein G(olf) alpha;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	GNAL

Image



Western blot analysis of extracts from LOVO cells, using GNAL antibody.



Immunohistochemistry analysis of paraffin-embedded human brain tissue using GNAL antibody.