




GNAT1 Antibody

Product Code	CSB-PA791645
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P11488
Immunogen	Synthesized peptide derived from internal of Human GNAT1.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Specificity	The antibody detects endogenous levels of total GNAT1 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	GBT1; GNATR; guanine nucleotide-binding protein G; transducin alpha-1 ch;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	GNAT1

Image



Western blot analysis of extracts from COLO cells, using GNAT1 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.