




GNAT2 Antibody

Product Code	CSB-PA237874
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P19087
Immunogen	Synthesized peptide derived from internal of Human GNAT2.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous levels of total GNAT2 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	ACHM4; GNATC; guanine nucleotide binding protein (G protein); alpha transducing activity 2; guanine nucleotide-binding protein G(t) alpha-2
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	GNAT2

Image



Western blot analysis of extracts from COLO cells, using GNAT2 antibody.