




GNL1 Antibody

Product Code	CSB-PA295891
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P36915
Immunogen	Synthesized peptide derived from internal of Human GNL1.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total GNL1 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	GTP-binding protein HSR1; guanine nucleotide binding protein-like 1; guanine nucleotide-binding protein-like 1; HSR1; HSR1 GTP-binding protein
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	GNL1

Image



Western blot analysis of extracts from Jurkat cells, LOVO cells and K562 cells, using GNL1 antibody.