




GNL3L Antibody

Product Code	CSB-PA106099
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NVN8
Immunogen	Synthesized peptide derived from internal of Human GNL3L.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total GNL3L protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	guanine nucleotide binding protein-like 3 (nucleolar)-like; guanine nucleotide-binding protein-like 3-like protein;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	GNL3L

Image



Western blot analysis of extracts from HeLa cells, 293 cells and HT-29 cells, using GNL3L antibody.