




HLA-DPA1 Antibody

Product Code	CSB-PA177845
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P20036
Immunogen	Synthesized peptide derived from N-terminal of Human HA2Q.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total HA2Q protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	DP(W3); DP(W4); DPA1; HLA class II histocompatibility antigen; DP alpha 1 chain
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	HLA-DPA1

Image



Western blot analysis of extracts from Jurkat cells, using HA2Q antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.