





HOXA11/HOXD11 Antibody

Product Code	CSB-PA913554
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P31270
Immunogen	Synthesized peptide derived from internal of Human HOXA11/D11.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous levels of total HOXA11/D11 protein.
Tested Applications	ELISA,WB,IF;WB:1:500-1:3000,IF:1:100-1:500
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Homeobox protein Hox-A11; Hox-1l; HOXA11; HOX1l;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	HOXA11/HOXD11

Image



Western blot analysis of extracts from RAW264.7 cells, using HOXA11/D11 antibody.



Immunofluorescence analysis of HUVEC cells,
using HOXA11/D11 antibody.