




HOXA5 Antibody

Product Code	CSB-PA184271
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P20719
Immunogen	Synthesized peptide derived from N-terminal of Human HXA5.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total HXA5 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	homeo box 1C; homeo box A5; homeobox A5; Homeobox protein Hox-1C; Homeobox protein Hox-A5
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	HOXA5

Image



Western blot analysis of extracts from Jurkat cells, using HXA5 antibody.