




HOXA5 Antibody

Product Code	CSB-PA184271
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P20719
Immunogen	Synthesized peptide derived from N-terminal of Human HXA5.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total HXA5 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	homeo box 1C; homeo box A5; homeobox A5; Homeobox protein Hox-1C; Homeobox protein Hox-A5
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	HOXA5

Image



Western blot analysis of extracts from Jurkat cells, using HXA5 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.