




MSX2 Antibody

Product Code	CSB-PA030362
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P35548
Immunogen	Synthesized peptide derived from Internal of Human MSX2.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Specificity	The antibody detects endogenous levels of total MSX2 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	CRS2; FPP; Homeobox protein Hox-8; Homeobox protein MSX-2; HOX8
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MSX2

Image



Western blot analysis of extracts from COLO205 cells, using MSX2 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.