




IQCB1 Antibody

Product Code	CSB-PA997466
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15051
Immunogen	Synthesized peptide derived from internal of Human IQCB1.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous levels of total IQCB1 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	IQ motif containing B1; NPHP5; PIQ; SLSN5;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	IQCB1

Image



Western blot analysis of extracts from K562 cells, using IQCB1 antibody.