




AASDHPPT Antibody

Product Code	CSB-PA299054
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NRN7
Immunogen	Synthesized peptide derived from internal of Human AASDHPPT.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total AASDHPPT protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	4-phosphopantetheinyl transferase; AASD-PPT; ADPPT; alpha-aminoadipic semialdehyde dehydrogenase-phosphopantetheinyl transferase; EC 2.7.8.-
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	AASDHPPT

Image



Western blot analysis of extracts from COLO cells, using AASDHPPT antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.