





DCT Antibody

Product Code	CSB-PA934850
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P40126
Immunogen	Synthesized peptide derived from internal of Human DCT.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous levels of total DCT protein.
Tested Applications	ELISA,WB,IF;WB:1:500-1:3000,IF:1:100-1:500
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	L-dopachrome tautomerase [Precursor]; DCT; DT; EC 5.3.3.12; L-dopachrome Delta-isomerase
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	DCT

Image



Western blot analysis of extracts from K562 cells, using DCT antibody.



Immunofluorescence analysis of A549 cells, using DCT antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.