




LILRA2 Antibody

Product Code	CSB-PA940104
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N149
Immunogen	Synthesized peptide derived from internal of Human LILRA2.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total LILRA2 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Leukocyte immunoglobulin-like receptor subfamily A member 2 [Precursor]; Leukocyte immunoglobulin-like receptor 7; LIR-7; Immunoglobulin-like transcript 1; ILT-1
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	LILRA2

Image



Western blot analysis of extracts from HepG2 cells, using LILRA2 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.