




LILRB4 Antibody

Product Code	CSB-PA297929
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NHJ6
Immunogen	Synthesized peptide derived from internal of Human LILRB4.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total LILRB4 protein.
Tested Applications	ELISA,IF;IF:1:100-1:500
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Leukocyte immunoglobulin-like receptor subfamily B member 4 ; Leukocyte immunoglobulin-like receptor 5; LIR-5; Immunoglobulin-like transcript 3; ILT-3
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	LILRB4

Image



Immunofluorescence analysis of A549 cells, using LILRB4 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.