



# SLC27A4 Antibody

<b>Product Code</b>	CSB-PA193776
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q6P1M0
<b>Immunogen</b>	Synthesized peptide derived from internal of Human SLC27A4.
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human, Mouse
<b>Specificity</b>	The antibody detects endogenous levels of total SLC27A4 protein.
<b>Tested Applications</b>	ELISA, IF; IF: 1:100-1:500
<b>Form</b>	Rabbit IgG in phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
<b>Purification Method</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Long-chain fatty acid transport protein 4; EC 6.2.1.-; Fatty acid transport protein 4; FATP-4; Solute carrier family 27 member 4
<b>Product Type</b>	Polyclonal Antibody
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	SLC27A4

## Image



Immunofluorescence analysis of A549 cells, using SLC27A4 antibody.

**Usage** For Research Use Only. Not for use in diagnostic or therapeutic procedures.