




SLC27A4 Antibody

Product Code	CSB-PA193776
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6P1M0
Immunogen	Synthesized peptide derived from internal of Human SLC27A4.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous levels of total SLC27A4 protein.
Tested Applications	ELISA,IF;IF:1:100-1:500
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Long-chain fatty acid transport protein 4; EC 6.2.1.-; Fatty acid transport protein 4; FATP-4; Solute carrier family 27 member 4
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SLC27A4

Image



Immunofluorescence analysis of A549 cells, using SLC27A4 antibody.