




LDLRAD1 Antibody

Product Code	CSB-PA951267
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5T700
Immunogen	Synthesized peptide derived from internal of Human LDLRAD1.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total LDLRAD1 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Low-density lipoprotein receptor class A domain-containing protein 1; LDLRAD1;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	LDLRAD1

Image



Western blot analysis of extracts from Jurkat cells and HepG2 cells, using LDLRAD1 antibody.