




LRP10 Antibody

Product Code	CSB-PA220915
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7Z4F1
Immunogen	Synthesized peptide derived from internal of Human LRP10.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous levels of total LRP10 protein.
Tested Applications	ELISA,IHC;IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Low-density lipoprotein receptor-related protein 10;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	LRP10

Image



Immunohistochemistry analysis of paraffin-embedded human brain tissue using LRP10 antibody.