




MED18 Antibody

Product Code	CSB-PA275975
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BUE0
Immunogen	Synthesized peptide derived from internal of Human MED18.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous levels of total MED18 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	FLJ20045; mediator 18; mediator of RNA polymerase II transcription 18; p28b;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MED18

Image



Western blot analysis of extracts from Jurkat cells and COLO cells, using MED18 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.