




MBTPS2 Antibody

Product Code	CSB-PA905280
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43462
Immunogen	Synthesized peptide derived from internal of Human MBTPS2.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Specificity	The antibody detects endogenous levels of total MBTPS2 protein.
Tested Applications	ELISA, IF; IF: 1:100-1:500
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Membrane-bound transcription factor site-2 protease; EC 3.4.24.85; S2P endopeptidase; Site-2 protease; Sterol regulatory element-binding proteins intramembrane protease
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MBTPS2

Image



Immunofluorescence analysis of A549 cells, using MBTPS2 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.