




MGST1 Antibody

Product Code	CSB-PA211410
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P10620
Immunogen	Synthesized peptide derived from internal of Human MGST1.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total MGST1 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	EC 2.5.1.18; Microsomal glutathione S-transferase 1; Microsomal GST-1; Microsomal GST-I;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MGST1

Image



Western blot analysis of extracts from HeLa cells, using MGST1 antibody.