




MRPS27 Antibody

Product Code	CSB-PA022313
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92552
Immunogen	Synthesized peptide derived from internal of Human MRPS27.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Specificity	The antibody detects endogenous levels of total MRPS27 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	mitochondrial 28S ribosomal protein S27; mitochondrial ribosomal protein S27; MRP-S27; RT27; S27mt
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MRPS27

Image



Western blot analysis of extracts from HeLa cells and HUVEC cells, using MRPS27 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.