




MRPL42 Antibody

Product Code	CSB-PA085001
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y6G3
Immunogen	Synthesized peptide derived from internal of Human MRPS32.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total MRPS32 protein.
Tested Applications	ELISA,IHC;IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Mitochondrial 28S ribosomal protein S32; mitochondrial; S32mt; MRP-S32;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MRPL42

Image



Immunohistochemistry analysis of paraffin-embedded human placenta tissue using MRPS32 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.