




MRPS34 Antibody

Product Code	CSB-PA944227
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P82930
Immunogen	Synthesized peptide derived from C-terminal of Human MRPS34.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Specificity	The antibody detects endogenous levels of total MRPS34 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	mitochondrial 28S ribosomal protein S34; mitochondrial ribosomal protein S34; MRP-S12; MRP-S34; RT34
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MRPS34

Image



Western blot analysis of extracts from HeLa cells, HUVEC cells, HT-29 cells, A549 cells and Jurkat cells, using MRPS34 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.